


Figure 4: Paired-End Sequencing and Alignment




Paired-end sequencing enables both ends of the DNA fragment to be sequenced. Because the distance between each paired read is known, alignment algorithms can use this information to map the reads over repetitive regions more precisely. This results in much better alignment of the reads, especially across difficult-to-sequence, repetitive regions of the genome.

End-to-End Solution

Only Illumina NGS provides a fully supported solution from DNA to results, with specialized sample prep choices for the application you are working on, to robust and proven sequencing reagents, and a wide range of simple data analysis tools (Figure 5).

Figure 5: Illumina’s End-to-End NGS Workflow on the MiSeq Benchtop Sequencer



NGS workflow on the MiSeq system includes sample preparation, massively parallel sequencing, automated data analysis, and cloud-enabled data analysis, storage, and sharing.

Take Your Research to the Next Level

The advent of NGS has enabled researchers to study biological systems at a level never before possible. With clear benefits over Sanger-based CE sequencing, next-generation sequencing can transform your microbiology research, opening new avenues to explore. To identify the sequencing platform that is optimal for your research needs, visit www.illumina.com.

From Innovation to Publication

As NGS technology continues to evolve, researchers are making fascinating discoveries in a number of biological fields, unlocking answers never before possible in all fields of research. As a result, there has been an explosion in the number of peer-reviewed scientific publications, including over 4,500 featuring Illumina sequencing technology. Selected recent examples relevant to microbiology are listed below.

Whole-Genome Sequencing

1. Toprak E, Veres A, Michel JB, Chait R, Hartl DL, et al. (2011) Evolutionary paths to antibiotic resistance under dynamically sustained drug selection. *Nat Genet* 44: 101–105.
2. Chua, KYL, Seemann T, Harrison PF, Monagle S, Korman TM, et al. (2011) The dominant Australian community-acquired methicillin-resistant *Staphylococcus aureus* clone ST93-IV [2B] is highly virulent and genetically distinct. *PLoS ONE* 6:
3. Loman NJ, Misra RV, Dallman TJ, Constantinidou C, Gharbia SE, et al. (2012) Performance comparison of benchtop high-throughput sequencing platforms. *Nat Biotechnol* 30(5): 434–9.

De novo Sequencing

4. Chitsaz H, Yee-Greenbaum JL, Tesler G, Lombardo MJ, Dupont CL, et al. (2011) Efficient de novo assembly of single-cell bacterial genomes from short-read data sets. *Nat Biotechnol* 29: 915–921.
5. Rodrigue S, Malmstrom R, Berlin A, Birren B, Henn M, et al. (2009) Whole genome amplification and de novo assembly of single bacterial cells *PLoS One* (4)9 e6864.

Metagenomics

6. Caporaso JG, Lauber CL, Walkers, WA, Berg-Lyons D, Lozupone CA et al. (2011) Global patterns of 16S rRNA diversity at a depth of millions of sequences per sample. *Proc Natl Acad Sci USA* 108:4516–22.
7. Mackelprang, R, Waldrop MP, DeAngelis KM, David MM, Chavarria KL, et al. (2011) Metagenomic analysis of a permafrost microbial community reveals a rapid response to thaw. *Nature* 480: 368–371.

© 2013–2014 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAlIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NextSeq, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub No. 770-2012-023 Current as of 20 December 2013

FOR RESEARCH USE ONLY

